

Experience The Difference:

AMS SYSTEMS = LOWER COST OF OWNERSHIP



DURABLE REACTION CUVETTES


Significant Lower Cost per Test/Panel
Reduced Medical Waste

AMS DURABLE CUVETTES are molded from optical grade plastic material and are individually washed and tested before made available for the next test; a software driven QA/QC program tests each cuvette and alerts the operator to change the one that needs to be replaced. This program allows prolonged cuvette usage while maintaining the analyzer precision characteristics thus reducing tests cost and medical waste.



LOWER COST PER TEST

Only 200uL Reagent Volume Required
Smaller Reagent Containers



AMS Diagnostics reagents and assays are very competitively priced and packaged conveniently for small to medium users test volume. AMS Diagnostics proprietary reagent containers are designed to have less than 300 uL dead volume for minimal reagent waste. The LIASYS 450 requires only 200uL of reagent volume per test.



LESS DOWN-TIME

10 month MTBF
Extended Warranty Plan

MTBF (Mean Time Between Failures) measures down-time due to component(s) failures. AMS LIASYS 330 series system boast 10 month MTBF time, over 1,267 installations over 3 years time. AMS Diagnostics Extended Warranty Plan provides for scheduled preventative maintenance visits, virtually eliminating service urgencies. Our Service Engineers are available directly at 1 (866) 369-2553.

REQUEST A COMPARATIVE ROI ANALYSIS AND
START REALIZING SAVINGS IN YOUR LABORATORY BY ADOPTING THE

AMS LIASYS CLINICAL CHEMISTRY SERIES

1(866) 419-7839 (Toll Free)
Service 1 (866) 369-2553 (Toll Free)